


EYECRYL™ PLUS

NATURAL YELLOW ASPHERIC HYDROPHOBIC ACRYLIC
FOLDABLE IOL WITH DELIVERY SYSTEM



FEATURES

- OPTIMALLY HYDRATED TO ENSURE "NO GLISTENINGS"
- NATURAL CHROMOPHORE
- 360 DEGREE ADVANCED SQUARE EDGE
- NEGATIVE S.A. TO COMPENSATE POSITIVE S.A. OF CORNEA
- EXCELLENT FOLDING & UNFOLDING PROPERTIES
- CONSISTENT BEHAVIOR IN ALL ENVIRONMENTAL CONDITIONS



biotech

VISION CARE

SPECIFICATIONS

Model: ASHFY600

MATERIAL	Hydrophobic Acrylic containing Natural Chromophore									
OPTIC TYPE	Single Piece, 360° Square Edge with Aspheric Optic									
OPTIC SIZE	6.00 mm									
OVERALL SIZE	13.00 mm									
ANGULATION	0°									
ACD	5.28									
REFRACTIVE INDEX	1.48									
ESTIMATED A-CONSTANT	118.5									
RECOMMENDED OPTICAL A-CONSTANTS	<table border="1"> <tr> <td>SRK - T 118.6</td> <td>SRK - II 118.9</td> <td>Holl 1 Const SF : 1.63</td> </tr> <tr> <td colspan="2">HOFFER Q ACD : 5.43</td> <td>HAIGIS $\alpha 0:1.20, \alpha 1:0.40, \alpha 2:1.10$</td> </tr> <tr> <td>Barrett 1.67</td> <td>Hill-RBF 118.60</td> <td>Olsen ACD 4.47</td> </tr> </table>	SRK - T 118.6	SRK - II 118.9	Holl 1 Const SF : 1.63	HOFFER Q ACD : 5.43		HAIGIS $\alpha 0:1.20, \alpha 1:0.40, \alpha 2:1.10$	Barrett 1.67	Hill-RBF 118.60	Olsen ACD 4.47
	SRK - T 118.6	SRK - II 118.9	Holl 1 Const SF : 1.63							
	HOFFER Q ACD : 5.43		HAIGIS $\alpha 0:1.20, \alpha 1:0.40, \alpha 2:1.10$							
Barrett 1.67	Hill-RBF 118.60	Olsen ACD 4.47								
DIOPTER RANGE	+5.0 to +30.0 D (with 0.5 D steps)									
IMPLANTATION SITE	Capsular Bag									
STERILIZATION	Irradiation									
SHELF LIFE	4 years from date of manufacture									

