

CE
2460

biotech

VISION CARE

CATARACT SURGERY SOLUTIONS

© 2014 Biotech Vision Care, All rights reserved.

Eye Visc

Hydroxy Propyl Methyl Cellulose (HPMC) Ophthalmic Solution USP 2.0% W/V


- Highly retentive, Dispersive viscoelastic
- Unique ability to offer good space maintenance and excellent tissue protection throughout the procedure
- Possesses a high dynamic viscoelasticity and coating ability
- Improves the depth of anterior chamber allowing easy and safe implantation of IOLs
- Excellent ease of use
- Viscosity 4500-9500 mPas

Indications:

- Eye Visc is indicated when ever protection & lubrication on the delicate cells & tissues are needed especially in ophthalmic procedure including Anterior Segment surgery, Cataract surgery, IOL implantation, Corneal transplantation & Glaucoma surgery.

Packaging:

- Eyevisc is supplied 2.0 ml in sterile PFS



biotech

VISIONCARE

TECNOIMAGEN

www.tecnoimagen.com.ar | 0810 333 8273

BIOTECH HEALTHCARE S.A. | Calle 12 de Abril 1000, Ciudad de Buenos Aires, Argentina | 1090
Calle 12 de Abril 1000, Ciudad de Buenos Aires | 1090
Email: info@biotechhealthcare.com | info@biotechhealthcare.com | www.biotechhealthcare.com